

Wondfo®

Blood Gas Analyzer

Test Items

Basic blood-gas:	pH, pCO ₂ , pO ₂
Haematology:	Hct
Electrolyte:	K ⁺ , Na ⁺ , Ca ⁺⁺ , Cl ⁻

Test Packages

Test Card Parameters	BGDR-3	BGDR-8	BGDR7-Ca	BGDR7-Cl	BGDR6	BGDR3	BGDR4
pH	*	*	*	*	*		
pCO ₂	*	*	*	*	*		
pO ₂	*	*	*	*	*		
K ⁺		*	*	*	*	*	*
Na ⁺		*	*	*	*	*	*
Ca ⁺⁺		*		*			*
Cl ⁻		*	*		*	*	*
Hct		*	*	*			

Wondfo®

Guangzhou Wondfo Biotech Co., Ltd.

Add: No. 8 Lizhishan Road, Science City, Luogang District, 510663, Guangzhou, P.R. China
Tel: +86-20-32296083 Fax: +86-20-32296063
E-mail: sales@wondfo.com.cn Websites: www.wondfo.com.cn

Version 2017/08/22

Wondfo®



BGA-101

Blood Gas Analyzer


- Quickly
- Portable
- Accurate
- Efficient

GUANGZHOU WONDFO BIOTECH CO., LTD.


Blood Gas Analyzer



Operation



Contents



Test Card

- Built-in injection slot, no cross contamination
- Dry chemical method, no clogging risk
- Less blood, suitable for newborns
- Just need 30s for detection
- Arterial blood, venous blood, capillary blood
- Good correlation to large laboratories

